medi Soft OA

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medi Soft OA


Donning instructions

Knee brace for strain relief and guidance


Note: It is best if you sit on the edge of a chair or a stool when putting on the brace.

The relief mechanism should always be in the starting position, even when it is not being worn. Important: The button for this is located at the white spot.

Open all quick closures on the front of the brace by lightly turning them. (Image 1).

Now open the Velcro closures on the brace.

Bend the affected leg to an approx. 70 degree angle and slip the brace over it (image 2).

Pull the brace upwards until the fabric touches the lower edge of the kneecap, making sure that the kneecap is positioned in the centre of the opening. The central hinge of the lateral hinge splints should now be approximately at the level of the upper edge of the kneecap (image 3).


Now close the upper Velcro closures. Then close the quick closures from strap 1 to strap 4 (image 4).

Push or pull the relief mechanism along the scale, either upwards or downwards. Each level is confirmed with a snap sound. The highest level of relief is attained at level 9.

Note: The pre-installation is intended as a basis for orientation and can be individually adjusted. Walk a few steps to sure that the brace is fitted correctly to your leg. Adjust the fit if necessary (image 5).

In order to reduce the pressure, unlock the mechanism and push it down the scale. Then activate it again (image 6).

medi tip

At www.medi.biz/oa you can find many helpful videos on how to fit braces correctly. There are also step-by-step instructions here.

Intended purpose: medi Soft OA is a knee brace for relief and guidance of the knee joint

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medi GmbH & Co. KG · Medicusstraße 1 · 95448 Bayreuth · Germany · www.medi.de/en