
TWO-LAYER COMPRESSION KIT Instructions For Use

(ENG)	Two-Layer Compression
(AR)	ضغط ثنائي الطبقات
(CS)	Dvouvrstvá komprese
(DA)	To-lags kompres
(DE)	2-Schicht-Komprimierung
(ES)	Compresión de dos capas
(FI)	2-kerroksinen kompressiosidos
(FR)	Compression double couche
(HE)	קומפרסיה דו-שכבתית
(HR)	Dvoslojna kompresa
(IT)	Compressione a due strati
(MS)	Mampatan Dua Lapisan
(NL)	Tweelaagse compressie
(NO)	To-lags kompresjon
(PL)	Dwuwarstwowy system uciskowy
(PT)	Compressão de duas camadas
(RU)	Двухслойная компрессия
(SK)	Dvojvrstvá kompresia
(SL)	Dvoslojna kompresija
(SV)	Tvålagers kompressionskit
(ZH-CN)	两层压缩
(ZH-TW)	兩層壓縮

Description

These instructions are for CoFlex TLC, CoFlex TLC XL and CoFlex TLC Lite.

CoFlex TLC (*Standard & Lite*) is a two-step non-latex compression system. Layer 1 is a foam roll, with odor and itch control properties, that provides comfort and absorbs drainage. Layer 2 is a short stretch compression bandage that provides the proper levels of sustained working pressures to promote venous return. This system is designed to work together and provide continuous restorative compression for up to one week. A nylon stocking is included for ease of use and patient comfort.

Indications for Use

CoFlex TLC (*Standard & Lite*) is designed to treat venous leg ulcers and other related conditions.

CoFlex TLC (*Standard*) can be used for patients with an ABI greater than or equal to 0.8. CoFlex TLC (*Lite*) can be used for patients with an ABI greater than or equal to 0.5.

Contraindications

CoFlex TLC (*Standard*): ABI less than 0.8.
CoFlex TLC (*Lite*): ABI less than 0.5.

Precautions & Warnings

- Should be used under the care of a licensed healthcare professional.
- Measure the Ankle Brachial Index prior to application.
- CoFlex TLC (*Standard & Lite*) Layer 1 is not a primary wound dressing. Follow appropriate protocol for the wound condition before applying.
- CoFlex TLC (*Standard & Lite*) Layer 2 may initially be applied with less tension allowing the patient to adjust to or tolerate the short stretch compression system.
- At the discretion of a healthcare professional, patients or their care providers can be trained to apply the bandage.
- Instruct patient to contact physician or nurse if they experience pain, numbness, tingling, discoloration or if swelling occurs.
- Product is not made with natural rubber latex.
- Single use item.
- If redness or irritation occurs, discontinue use.

Directions for Use

General Instructions:
Position the foot at a 90° angle.

Layer 1 - Foam Layer:

- (1) Begin wrapping Layer 1 (with colored side against the skin) on the top of the foot, at the base of the fifth digit in a spiral technique (figure 1).
- (2) Continue spiral wrapping across the top of the foot then proceed to wrap 1 figure of 8 around the back of the heel/ankle. The plantar surface of the heel may be partially exposed (figure 2).
- (3) Continue spiral wrapping up the leg. Apply with a 50% overlap using just enough tension to conform to the shape of the leg, ending just below the knee (figure 3).
- (4) Cut or hand-tear excess material and tape to keep the bandage in place (figure 4).

Layer 2 - Compression Layer:

- (1) Begin wrapping Layer 2 on the top of the foot, at the base of the fifth digit in a spiral technique (figure 5).
- (2) Hold the bandage with thumb and stretch until the ovals become circles. Begin spiral wrapping across the top of the foot, then proceed to wrap two figure of 8's around the back of the heel/ankle, completely covering the heel and foot. Continue to maintain even tension by stretching the visual indicators from ovals to circles (figure 6). For 2" size, indicators are not included - stretch to the appropriate level.

- (3) Continue to spiral wrap up the leg, covering Layer 1 with 50% overlap while maintaining even tension (figure 7).
- (4) Cut or hand-tear excess material (figure 8).
- (5) Gently press over entire application to ensure a strong bond & reduce slippage (figure 9).
- (6) Apply the nylon stocking over the system for ease of movement and patient comfort in footwear, under clothes and on bed sheets (figure 10).

Removal

CoFlex TLC (*Standard & Lite*) may be removed by unwrapping or with bandage scissors.

Storage

For best results, store CoFlex TLC (*Standard & Lite*) in its original packaging at room temperature until use. Avoid excessive heat and humidity. Refer to the label on each box for expiration date.

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CoFlex is a registered trademark of Andover Healthcare, Inc., a subsidiary of Milliken Healthcare Products, LLC.
Advena, Ltd.
Andover Healthcare, Inc.