

Art. nr.	Navn	Mål i cm							
		Antall celler	A Lår	B Over kne	C Legg	D Ankel	E Vrist	F Forfot	L Lengde
29F0836	Lympha Press 8 3-65	8	78	71	65	50	54	38	65
29F0837	Lympha Press 8 3-75	8	81	74	66	50	50	38	75
29F0838	Lympha Press 8 3-85	8	86	80	68	53	51	38	85
29F0839	Lympha Press 4 2-50	4	61	54	46	43	48	26	50
29F0840	Lympha Press 4 3-65	4	78	71	65	50	54	38	65
29F0841	Lympha Press 4 3-75	4	81	74	66	50	50	38	75
29F0842	Lympha Press 4 3-85	4	86	80	68	53	51	38	85
29F0843	Lympha Press 4 4-65	4	105	98	86	67	60	41	65
29F0844	Lympha Press 4 4-75	4	107	99	86	72	60	41	75
29F0845	Lympha Press 4 4-85	4	109	105	93	67	60	43	85
29F0846	Personal leg 2-50	4	61	54	46	43	48	26	50
29F0847	Personal leg 2-65	4	63	55	46	43	48	26	65
29F0848	Personal leg 2-75	4	66	56	46	43	48	26	75
29F0849	Personal leg 2-85	4	70	60	47	43	48	26	85
29F0850	Personal leg 3-75	4	81	74	66	50	50	38	75