



Art. nr.	Navn	Omkrets i cm				Mansjett lengde	Fotlengde	Antall kammer
		Øverst på mansjett		Nede ved ankel				
		min	max	min	max			
		A	B	C	D			
29F0059	FL benmansjett 2-75	60	70	40	50	75	30	12
29F0061	FL benmansjett 2-85	66	76	40	50	85		
29F0803	FL benmansjett 2-95	70	80	40	50	95		
29F0804	FL benmansjett 2-105	74	84	40	50	105		
29F0805	FL benmansjett 3-75	66	76	43	53	75	34	
29F0060	FL benmansjett 3-85	70	80	43	53	85		
29F0062	FL benmansjett 3-95	74	85	43	53	95		
29F0083	FL benmansjett 4-75	78	90	56	68	75	36	
29F0806	FL benmansjett 4-85	88	96	56	68	85		
29F0807	FL benmansjett 4-95	90	108	56	68	95		