


Prescription pad medi orthopaedics

Knee braces

M.4 X-lock	M.4 s	M.4s OA	Collamed
4-point hard-frame brace with extension / flexion limitation	Functional rigid brace with physioglide hinge for flexion/ extension limitation	4-point knee orthosis for control of anterior and posterior tibial translation. 3-point off-loader for varus or valgus pressure relief	Long, functional soft brace with flexion/ extension limitation, made of medi Airtex-material
 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>
Stabimed	medi PTS	medi PT control	protect.PT soft
Short, functional soft brace with flexion/ extension limitation, made of medi Airtex-material	Posterior - Tibial - Support, knee immobilisation support	Soft brace for patella tracking control with flexion/extension limitation, made of medi Airtex-material	Soft brace for patellar tracking control
 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>	 <input type="checkbox"/>

Measurement



circumference of thigh 15 cm above the centre of the patella (measured at the medial side of the thigh)

M.4 X-lock, M.4s, M.4s OA, Collamed, Stabimed, medi PT control, protect.PT soft


XS	31 – 37 cm	circ.
S	37 – 45 cm	circ.
M	45 – 51 cm	circ.
L	51 – 57 cm	circ.
XL	57 – 65 cm	circ.
XXL	65 – 74 cm	circ.

medi PTS

45 cm	length
55 cm	length
65 cm	length


Supports

Achimed	Levamed	Epicomed	Genumedi
Achilles tendon support with integrated anatomically formed silicone pad	Ankle support with silicone cushions	Elbow soft support with silicone support pads and tension strap	Knee support with patella silicone ring



Information

To order by Credit Card:
phone 01432 373500



Genumedi

(f) measurement: circumference of thigh 15 cm (6 inches) above centre of patella

(d) measurement: circumference 5 cm (2 inches) below centre of patella

		I	II	III	IV	V	VI	VII
standard (d)	circumference	25 - 28 cm	28 - 31 cm	31 - 34 cm	34 - 37 cm	37 - 40 cm	40 - 43 cm	43 - 46 cm
standard (f)	circumference	37 - 40 cm	40 - 43 cm	43 - 46 cm	46 - 49 cm	49 - 52 cm	52 - 55 cm	55 - 58 cm
extra wide (f)	circumference				49 - 52 cm	52 - 55 cm	55 - 59 cm	58 - 61 cm



Achimed / Levamed

size	I	II	III	IV	V	VI
circumference	17 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 26 cm	26 - 28 cm	28 - 30 cm

measurement:
circumference of ankle
at narrowest point
above ankle bone



Epicomed

size	1	2	3	4	5	6
circumference	17-20 cm	20-23 cm	23-27 cm	27-30 cm	30-32 cm	32-35 cm

measurement:
circumference 5 cm
below centre of elbow
joint

