


# Information och måttabell

## Rund- & flatstickat, arm standardstorlekar



### mediven harmony

Armstrumpor omkrets (cm)							
Storlek	I	II	III	IV	V	VI	VII
cG x-vid	28,5-31,0	30,5-33,5	33,0-36,0	35,0-38,5	37,0-40,5	39,5-43,0	41,5-45,5
cG	25,5-28,0	28,0-30,5	30,0-33,0	32,0-35,0	33,5-36,5	36,0-39,5	38,0-41,5
cF x-vid	25,5-28,0	27,5-30,0	30,0-33,0	32,0-35,0	34,0-37,0	36,0-39,5	38,5-42,0
cF	23,5-25,5	25,0-27,5	27,5-30,0	29,5-32,0	31,0-34,0	33,0-36,0	35,0-38,5
cE x-vid	24,0-26,0	26,5-29,0	28,5-31,0	30,0-33,0	31,5-34,5	33,5-36,5	35,5-39,0
cE	23,0-25,0	24,5-26,5	26,0-28,5	27,5-30,0	29,0-31,5	30,5-33,5	32,5-35,5
cD x-vid	24,0-26,0	25,0-27,5	26,0-28,5	28,0-30,5	29,5-32,0	31,0-34,0	33,5-36,5
cD	22,0-24,0	23,0-25,0	24,0-26,0	25,0-27,5	27,0-29,5	28,5-31,0	30,5-33,5
cC <sub>1</sub> = 100%*	15,0-16,5	16,5-18,0	18,0-19,5	19,5-21,0	21,0-23,0	23,0-25,0	25,0-27,5
cC	13,5-14,5	14,5-15,5	16,0-17,5	17,5-19,0	19,0-20,5	20,5-22,5	22,5-24,5
cA	15,0-18,0	16,0-19,0	17,0-20,0	19,0-22,0	20,0-23,0	22,0-25,0	24,0-27,0

Längd (cm)	
Standardlängd	41 – 47

\* Mätpunkt cC<sub>1</sub> är 5 cm ovanför cC

### Handskar omkrets (cm)

Storlek	II	III	IV	V	VI	VII
cC <sub>1</sub> = 100%*	15,0-18,0	17-20	19-21,5	21-23,5	23-26	25,5-30
cC	13,5-15,5	15-18	17-20	19-21	20-23,5	22,5-27,5
cB	15,5-20	17,5-21	19,5-23	20,5-24	22,5-26	24,5-30
cA	15-19	17-20	19-20	20-23	22-25	24-29

### Fingrar omkrets (cm)

Storlek	II	III	IV	V	VI	VII	
Tumme	Bas	7,0	7,5	8,0	8,5	9,0	10,5
	Topp	6,0	6,5	6,5	7,0	7,5	8,0
Pekfinger	Bas	6,5	7,0	7,5	8,0	8,5	9,5
	Topp	5,0	5,5	6,0	6,0	6,5	7,0
Långfinger	Bas	6,5	7,0	7,5	8,0	8,5	9,5
	Topp	5,0	5,5	6,0	6,0	6,5	7,0
Ringfinger	Bas	6,0	6,5	7,0	7,0	7,5	9,0
	Topp	5,0	5,0	5,5	5,5	6,0	6,5
Lillfinger	Bas	5,5	6,0	6,0	6,5	7,0	8,0
	Topp	4,5	4,5	5,0	5,5	5,5	6,0

### Handskar längd (cm)

AB	3,5-4,5
AC	10-12
AC <sub>1</sub>	15-16

### Fingrar längd (cm)

Tumme	3,5-4,5
Pekfinger	4,5-5,5
Långfinger	5,5-6,5
Ringfinger	4,5-5,5
Lillfinger	3,5-4,5

\* Mätpunkt cC<sub>1</sub> är 5 cm ovanför cC

### mediven esprit

Armstrumpor omkrets (cm)					
Storlek	I	II	III	IV	V
cA	17-19	18-20	19-22	21-24	23-27
cC	15-17	17-19	19-21	21-23	23-25
cE	25-27,5	27-29,5	29-32	31-34	33-37
cG	28-32	30-34	32-37	35-40	38-44

Standardlängd cG: 40-44cm, lång: 44-48cm

Handskar omkrets (cm)					
Storlek	I	II	III	IV	V
A	17,5-19,5	18,5-20	20-22	22-24	24-27
B	18-19	19-21,5	21,5-23,5	23,5-25,5	25,5-28,5
C	15-17	17-19	19-21,5	21-23	23-25

### Handskar längd (cm)

AB	3,5-4,5
AC	10-12
AC <sub>1</sub>	15-16

### Fingrar längd (cm)

Tumme	3,5-4,5
Pekfinger	4,5-5,5
Långfinger	5,5-6,5
Ringfinger	4,5-5,5
Lillfinger	3,5-4,5