

Patellar braces



medi Patella tendon support
Patella tendon support

Genumedi®
Knee support with patella silicone ring

Genumedi® PT
Knee support promotes muscular balance in cases of patellar lateralisation

Genumedi® PSS
Knee support for relieving strain on the patellar tendon entheses

protect.PT soft
Soft brace for patellar tracking control

medi PT control
Soft brace for patella tracking control with flexion / extension limitation

- Pressure on the quadriceps tendon improves proprioceptive performance and thus reinforces muscular support of the knee joint

- The massage effect achieved by the compressive fabric and silicone insert improves the blood circulation and aids absorption of haematomas and oedemas
- The seamless, tailored knee support improves proprioception and improves muscular support (stabilisation of the knee) Joint effusion and swelling Irritative conditions

also available: extra wide and sand

- Positive effect on the muscular structures that guide the patella to compensate for imbalance
- Reduces muscle tone of vastus lateralis
- Activates vastus medialis
- Laterally reinforced silicone ring gently guides the patella for reliable movement in the femoral groove
- Proprioceptive compression for reliable guidance and stabilisation of the joints

- Two-part pad system for easy patellar guidance and specific reduction of pressure and tension at the tip of the patella and the tibial tuberosity
- Offset strap system for balanced pressure distribution
- Pressure exerted by the integrated silicone dot coating stimulates proprioception to increase quadriceps muscle tone
- Proprioceptive compression for reliable guidance and stabilisation of the joints

- Correction of lateral patella maltracking due to a visco-elastic stabilisation pad
- Soft brace with compressive and proprioceptive effect

- Collateral stabilisation of the knee joint
- Prevention of hyperextension and limitation of range of motion
- Correction of lateral patella mal-tracking
- Compressive and proprioceptive effect
- Patella stabilisation buttress that prevents lateral movement throughout full flexion and extension

size	code no.	size	code no.	size	code no. left	code no. right	size	code no.	size	code no. left	code no. right	size	code no. left	code no. right
universal	877	–	–	–	–	–	0	620 / 0	–	–	–	–	–	–
		I	613 / I	I	619 / I	618 / I	I	620 / I	1	P775 / 1	P774 / 1	1	889 / 1	888 / 1
		II	613 / II	II	619 / II	618 / II	II	620 / II	2	P775 / 2	P774 / 2	2	889 / 2	888 / 2
		III	613 / III	III	619 / III	618 / III	III	620 / III	3	P775 / 3	P774 / 3	3	889 / 3	888 / 3
		IV	613 / IV	IV	619 / IV	618 / IV	IV	620 / IV	4	P775 / 4	P774 / 4	4	889 / 4	888 / 4
		V	613 / V	V	619 / V	618 / V	V	620 / V	5	P775 / 5	P774 / 5	5	889 / 5	888 / 5
		VI	613 / VI	VI	619 / VI	618 / VI	VI	620 / VI	6	P775 / 6	P774 / 6	6	889 / 6	888 / 6
		VII	613 / VII	VII	619 / VII	618 / VII	VII	620 / VII						

