

# The OA Knee Brace Portfolio



OA grades

**protect.OA soft**  
soft brace for  
Osteoarthritis

**medi Soft OA light**  
soft knee orthosis for pressure  
relief and guidance

**medi Soft OA**  
soft knee orthosis for pressure  
relief and guidance

**M.3 s OA**  
single hinged functional  
knee orthosis for the relief  
of unicompartmental  
osteoarthritis of the knee

**The Game Changer**  
universal knee brace

**protect.4 OA**  
off loading knee orthosis  
for control of anterior and  
posterior tibial translation and  
varus or valgus pressure relief

**M.4® s OA comfort**  
4-point knee orthosis for control of  
anterior and posterior tibial  
translation. 3-point off-loader for  
varus or valgus pressure relief

Soft support for mild unicompartmental osteoarthritis. Tensioning straps provide varus / valgus relief.

Relief of the lateral or medial compartments through three-point principle

Relief of the lateral or medial compartments through three-point principle

Single hinge applies pressure to the medial or lateral aspect to correct varus / valgus alignment. Ideal for use where bilateral hinges cannot be tolerated.

The simple design makes it easy to use and helps to restore full range of motion from OA of the knee.

Lightweight alloy brace, provides more protection and control than the soft braces. Varus / valgus pressure applied through medial / lateral hinge.

Strong alloy brace that accurately mimics the glide / roll action of the knee joint whilst applying targeted varus / valgus pressure to relieve the affected compartment. Designed for long term use.

size	code no.	size	code no.	size	code no.	size	code no. left	code no. right	code no. left	code no. right	size	code no. left	code no. right	size	code no. left	code no. right
left varus / right valgus		left varus / right valgus		left varus / right valgus		varus	left	right	left	right	varus	left	right	varus	left	right
XS	P773 / 1	XS	915 / 1	XS	913 / 1	XS	909 / 1	908 / 1	917	916	XS	P051 / 1	P050 / 1	XS	905 / 1	904 / 1
S	P773 / 2	S	915 / 2	S	913 / 2	S	909 / 2	908 / 2			S	P051 / 2	P050 / 2	S	905 / 2	904 / 2
M	P773 / 3	M	915 / 3	M	913 / 3	M	909 / 3	908 / 3			M	P051 / 3	P050 / 3	M	905 / 3	904 / 3
L	P773 / 4	L	915 / 4	L	913 / 4	L	909 / 4	908 / 4			L	P051 / 4	P050 / 4	L	905 / 4	904 / 4
XL	P773 / 5	XL	915 / 5	XL	913 / 5	XL	909 / 5	908 / 5			XL	P051 / 5	P050 / 5	XL	905 / 5	904 / 5
XXL	P773 / 6	XXL	915 / 6	XXL	913 / 6	XXL	909 / 6	908 / 6			XXL	P051 / 6	P050 / 6	XXL	905 / 6	904 / 6
left valgus / right varus		left valgus / right varus		left valgus / right varus		valgus	left	right			valgus	left	right	valgus	left	right
XS	P770 / 1	XS	914 / 1	XS	912 / 1	XS	911 / 1	910 / 1			XS	P053 / 1	P052 / 1	XS	907 / 1	906 / 1
S	P770 / 2	S	914 / 2	S	912 / 2	S	911 / 2	910 / 2			S	P053 / 2	P052 / 2	S	907 / 2	906 / 2
M	P770 / 3	M	914 / 3	M	912 / 3	M	911 / 3	910 / 3			M	P053 / 3	P052 / 3	M	907 / 3	906 / 3
L	P770 / 4	L	914 / 4	L	912 / 4	L	911 / 4	910 / 4			L	P053 / 4	P052 / 4	L	907 / 4	906 / 4
XL	P770 / 5	XL	914 / 5	XL	912 / 5	XL	911 / 5	910 / 5			XL	P053 / 5	P052 / 5	XL	907 / 5	906 / 5
XXL	P770 / 6	XXL	914 / 6	XXL	912 / 6	XXL	911 / 6	910 / 6			XXL	P053 / 6	P052 / 6	XXL	907 / 6	906 / 6

