

Knee ligament braces



Increasing support

Genumedi®
Knee support for soft tissue compression

Genumedi® pro
Knee support with plastic hinge and straps

Stabimed®
Short functional soft knee brace with tool free flexion / extension limitation

Collamed®
Long functional soft knee brace with tool free flexion / extension limitation and 4-point-principle

M.4s® comfort
4-point rigid brace with physioglide hinge for extension / flexion limitation

M.4 X-lock®
4-point rigid brace with physioglide hinge and extension lock

M.4s® PCL dynamic
4-point rigid frame brace with extension / flexion limitation for treating injuries of the posterior cruciate ligament

High quality knitted knee support with open patella ring and popliteal Comfort Zone area. Clima Comfort and Clima Fresh technologies enables the wearer to remain comfortable all day.

**also available:
extra wide and sand**

High quality knitted knee support with open patella ring and popliteal Comfort Zone area. Clima Comfort and Clima Fresh technologies enables the wearer to remain comfortable all day. Modern easyglide hinges stabilise and support the natural roll and glide motion of the knee joint.

medi Airtex+ spacer knit with Adaptive and Pure technologies ensures that maximum comfort is achieved and compliance remains optimal. Physioglide tool free hinge mimics the knees natural roll and glide for optimal fit during gait.

medi Airtex+ spacer knit with Adaptive and Pure technologies ensures that maximum comfort is achieved and compliance remains optimal. Physioglide tool free hinge mimics the knees natural roll and glide for optimal fit during gait.

Low profile functional brace for ligamentous instabilities. Manufactured from aircraft grade aluminium ensures it is as light as possible to aid comfort. Physioglide hinge mimics the knees natural roll and glide for optimal fit during gait.

Waterproof pads to facilitate wear during washing. Simple to use lock in full extension. Low profile functional brace for ligamentous instabilities. Manufactured from aircraft grade aluminium ensures it is as light as possible to aid comfort. Physioglide hinge mimics the knees natural roll and glide for optimal fit during gait.

The product allows the conservative or postoperative treatment of patients with a (partial) rupture of the posterior cruciate ligament. The M.4s PCL dynamic protection pad holds the tibia in the desired position. The integrated tensioning dial enables fast and precise adjustment.

size	code no.	size	code no.	size	code no.	size	code no.	size	code no. left	code no. right	size	code no. left	code no. right	size	code no. left	code no. right
I	613 / I	1	617 / 1	XS	826 / 1	XS	824 / 1	XS	801 / C / 1	800 / C / 1	XS	811 / 1	810 / 1	XS	813 / 1	812 / 1
II	613 / II	2	617 / 2	S	826 / 2	S	824 / 2	S	801 / C / 2	800 / C / 2	S	811 / 2	810 / 2	S	813 / 2	812 / 2
III	613 / III	3	617 / 3	M	826 / 3	M	824 / 3	M	801 / C / 3	800 / C / 3	M	811 / 3	810 / 3	M	813 / 3	812 / 3
IV	613 / IV	4	617 / 4	L	826 / 4	L	824 / 4	L	801 / C / 4	800 / C / 4	L	811 / 4	810 / 4	L	813 / 4	812 / 4
V	613 / V	5	617 / 5	XL	826 / 5	XL	824 / 5	XL	801 / C / 5	800 / C / 5	XL	811 / 5	810 / 5	XL	813 / 5	812 / 5
VI	613 / VI	6	617 / 6	XXL	826 / 6	XXL	824 / 6	XXL	801 / C / 6	800 / C / 6	XXL	811 / 6	810 / 6	XXL	813 / 6	812 / 6
VII	613 / VII	7	617 / 7													



4 046114 157907 009UK11/01.2019