

duomed soft® 2easy - 2 piece compression system



Comprising - 2 calf sleeves and 2 anklets
 or anklet pack (2 anklets)

All sizes below held in stock for next working day delivery.

Patient name: _____


Date: _____

Doctor/Clinician: _____

Clinic or Pharmacy Address:

(Please specify the address where the patient will collect their hosiery)

Request for duomed soft® 2easy on FP10/GP10 prescription



SIZE (circumference in cms)	SMALL (2)	MEDIUM (3)	LARGE (4)	X LARGE (5)	XX LARGE (6)
cC CALF	28—34	32—38	36—42	40—46	44—50
cB ANKLE	19—21	22—24	25—27	28—30	31—34

1 PATIENT MEASUREMENTS	LEFT LEG	RIGHT LEG	2	= Size requested
	cC Widest part of calf			
cB Just above ankle bone			3	= Number of pairs required

4 Choose your compression class:

BS ccl1—14-17 mmHg

BS ccl 2—18-24 mmHg

5 Choose your style:

Below knee (2 anklets & 2 calf sleeves)

Anklet (2 anklets only)

6 Choose your colour:

Sand

Black

Your patient should present this sheet to their GP to obtain product on prescription

duomed soft® 2easy - 2 piece compression system



Order Codes	BS ccl1 14-17mmHg		BS ccl2 18-24 mmHg	
Colour	Sand	Black	Sand	Black
duomed soft 2easy calf sleeve and ankle (pair)	DT2E1/size	DT2E1/B/size	DT2E2/size	DT2E2/B/size
Anklet (pair)	DT2EA1/size	DT2EA1/B/size	DT2EA2/size	DT2EA2/B/size

duomed soft 2easy (pair of calf sleeves and anklets)						
BS compression class 1 - 14-17mmHg						
	PRODUCT CODE	2 (small)	3 (medium)	4 (large)	5 (x large)	6 (xx large)
Sand	DT2E1/	417-1880	417-1898	417-1906	417-1914	417-1922
Black	DT2E1/B/	417-1930	417-1948	417-1955	417-1963	417-1971
BS compression class 2 - 18-24mmHg						
Sand	DT2E2/	417-2102	417-2110	417-2128	417-2136	417-2144
Black	DT2E2/B/	417-2151	417-2169	417-2177	417-2185	417-2193

duomed soft 2easy (pair of anklets only)						
BS compression class 1 - 14-17mmHg						
	PRODUCT CODE	2 (small)	3 (medium)	4 (large)	5 (x large)	6 (xx large)
Sand	DT2EA1/	417-2003	417-2029	417-2011	417-2037	417-2045
Black	DT2EA1/B/	417-2052	417-2060	417-2078	417-2086	417-2094
BS compression class 2 - 18-24mmHg						
Sand	DT2EA2/	417-2201	417-2219	417-2227	417-2235	417-2243
Black	DT2EA2/B/	417-2250	417-2268	417-2276	417-2284	417-2292

For further information, please call 01432 373 500 to speak directly with our Customer Services Team or email enquiries.clinical@mediuk.co.uk