

juxtafit®

Patient Referral Form FP10/GP10

lower legging & ankle foot wrap

Your patient should present this sheet to their GP to obtain product on Prescription.

juxtafit lower leg

Product length: petite 28cm or standard 36cm.
Measured from medial ankle bone to knee crease in standing position.

juxtafit ankle foot wrap


Closed heel, foot length from base of toes 16cm/
interlocking (open heel), foot length from base of toes 13-25cm.

From: _____

Patient's name: _____

Clinician: _____ Date: _____

Contact Tel: _____




juxtafit lower legging/lower legging extra wide calf

		2 (S)	3 (M)	3 (M-XW)	4 (L)	4 (L-XW)	5 (XL)	5 (XL-XW)	6 (XXL)
Circumferences in cm									
cC		30-38	36-43	43-50	41-48	48-56	41-51	51-61	53-64
cB1		20-28	25-33	33-40	30-38	38-46	33-41	41-51	43-53
cB		15-23	20-28	20-28	25-33	25-33	30-38	30-38	33-41
Petite length 28cm (11")	Product Code	DT950/2	DT950/3	DT950/3X	DT950/4	DT950/4X	DT950/5	DT950/5X	DT950/6
	PIP Code	348-0548	348-0555	368/0782	348-0563	368-0790	348-0571	368-0808	348-0589
Required ✓									
Standard length 36cm (14")	Product Code	DT951/2	DT951/3	DT951/3X	DT951/4	DT951/4X	DT951/5	DT951/5X	DT951/6
	PIP Code	348-0597	348-0605	368-0816	348-0613	368-0824	348-0621	368-0832	348-0639
Required ✓									

juxtafit ankle foot wrap, interlocking (open) & closed heel

		2 (S)	3 (M)	4 (L)
Circumferences in cm				
Length A-Y		>9	>9	>9
cH		19-24	24-29	29-34
cB		>36	>36	>41
Open heel	Product Code	DT960/2	DT960/3	DT960/4
	PIP Code	348-0647	348-0654	348-0662
Required ✓				
Closed heel	Product Code	DT961/2	DT961/3	DT961/4
	PIP Code	348-0670	348-0688	348-0696
Required ✓				




juxtafit®

Patient Referral Form FP10/GP10

juxtafit® upper leg and knee piece

Your patient should present this sheet to their GP to obtain product on Prescription.



juxtafit upper leg and knee piece

		1(XS)	2(S)	3(M)	4(L)	5(XL)
Circumferences in cm						
cG		46-62	54-70	60-80	69-90	77-101
cE1		37-48	42-55	50-65	58-76	66-87
cD		27-35	31-41	35-46	40-52	44-58
35cm length Right	Product Code PIP Code	DT954R/1 404-1216	DT954R/2 404-1224	DT954R/3 404-1232	DT954R/4 404-1240	DT954R/5 404-1257
Required ✓						
35cm length Left	Product Code PIP Code	DT954L/1 404-1133	DT954L/2 404-1117	DT954L/3 404-1158	DT954L/4 404-1174	DT954L/5 404-1190
Required ✓						
45cm length Right	Product Code PIP Code	DT955R/1 378-6548	DT955R/2 378-6555	DT955R/3 378-6563	DT955R/4 378-6571	DT955R/5 378-6589
Required ✓						
45cm length Left	Product Code PIP Code	DT955L/1 378-6498	DT955L/2 378-6506	DT955L/3 378-6514	DT955L/4 378-6522	DT955L/5 378-6530
Required ✓						
55cm length Right	Product Code PIP Code	DT956R/1 378-6654	DT956R/2 378-6662	DT956R/3 378-6670	DT956R/4 378-6688	DT956R/5 378-6696
Required ✓						
55cm length Left	Product Code PIP Code	DT956L/1 378-6597	DT956L/2 378-6605	DT956L/3 378-6613	DT956L/4 378-6621	DT956L/5 378-6639
Required ✓						