

juxtalite®

Patient Referral Form FP10/GP10

juxtalite is listed under Part IXA of the Drug Tariff, under 'Venous Ulcer Compression System'

From: _____	
Patient's name: _____	
Clinician: _____	Date: _____
Contact Tel: _____	

Your patient should present this sheet to their GP to obtain product on Prescription.

IB-D Length of garment following contour of leg

This garment is available in 2 lengths and 8 sizes.

To select the correct garment, measure length and circumferences as follows:

LENGTH: Length of garment from B (ankle) to D (2 fingers width below crease of the knee) following contour of leg. Select best fit.

28cm = short, 33cm = long.

CIRCUMFERENCES: Circumference at ankle (B) and widest calf (C).

Circumferences in cm								
	2 (S)	3 (M)	3 (M-XW)	4 (L)	4 (L-XW)	5 (XL)	5 (XL-XW)	6 (XXL)
cC	26-36	33-44	44-54	41-51	51-61	37-48	48-61	53-64
cB	19-22	22-27	22-27	25-32	25-32	30-37	30-37	35-42

juxtalite Codes									
Product Code & Length	Small	Medium	Medium X Wide	Large	Large X Wide	X Large	X Large X Wide	XX Large	
Short (28cm)	PIP Code	394-0996	394-1002	394-1044	394-1010	394-1051	394-1028	394-1069	394-1036
	Product Code	DT301/2	DT301/3	DT301X/3	DT301/4	DT301X/4	DT301/5	DT301X/5	DT301/6
	Quantity								
Long (33cm)	PIP Code	394-0913	394-0921	394-0962	394-0939	394-0970	394-0947	394-0988	394-0954
	Product Code	DT300/2	DT300/3	DT300X/3	DT300/4	DT300X/4	DT300/5	DT300X/5	DT300/6
	Quantity								

The juxtalite pack contains 1 x juxtalite, 1 x BPS card, 2 x black compressive foot juxta undersocks & 1 x instructions for use.

