

mediven[®] compression socks

Your patient should present this sheet to their GP to obtain product on Prescription.

From:	_____
Patient's name:	_____
Measured by:	Date: _____
Contact Tel:	_____

Selecting the correct garment for your patient.

mediven [®] for men	
Closed toe	
1. Choose your RAL Compression class:	
<input type="checkbox"/> Class 1 18 - 21mmHg	<input type="checkbox"/> Class 2 23 - 32mmHg

mediven [®] active	
Closed toe	
1. Choose your RAL Compression class:	
<input type="checkbox"/> Class 1 18 - 21mmHg	<input type="checkbox"/> Class 2 23 - 32mmHg

2. Choose your Style:			
Below knee (pair)			
Length	<input type="checkbox"/> Standard	<input type="checkbox"/> Petite	
Standard length ID (see reverse) - 39 - 44cm			
Petite length ID (see reverse) - 33 - 38cm			
Order codes - Class 1		Order codes - Class 2	
Standard	Petite	Standard	Petite
DT342	DT343	DT340	DT341

2. Choose your Style:			
Below knee (pair)			
Length	<input type="checkbox"/> Standard	<input type="checkbox"/> Petite	
Standard length ID (see reverse) - 43-49cm			
Petite length ID (see reverse) - 37-42cm (black only)			
Order codes - Class 1		Order codes - Class 2	
Standard	Petite	Standard	Petite
DT180	DT181	DT280	DT281

3. Choose your colour:			
<input type="checkbox"/> Black (code/5)	<input type="checkbox"/> Navy (code/3)	<input type="checkbox"/> Grey (code/2)	

3. Choose your colour:			
<input type="checkbox"/> Black (code/5)	<input type="checkbox"/> Navy (code/3)	<input type="checkbox"/> Grey (code/2)	
<input type="checkbox"/> White (code/7)	<input type="checkbox"/> Brown (code/6)		

4. Choose your size:							
<input type="checkbox"/> I	<input type="checkbox"/> II	<input type="checkbox"/> III	<input type="checkbox"/> IV	<input type="checkbox"/> V	<input type="checkbox"/> VI	<input type="checkbox"/> VII	
See over for measurement chart							

5. Pairs required: <input type="text"/>	6. DT code if known: <input type="text"/>
--	--

mediven[®] compression socks

TIPS FOR MEASURING:


- 1) Measure the circumference at the narrowest part of the ankle (B).
- 2) Measure the circumference at the widest part of the calf (C).
- 3) Measure the length of the limb at the back of the heel (flat on the floor) to a point 2 finger widths down from the knee flexure.

EXAMPLE:

Circumference B= 24cm, circumference C= 40cm, therefore could be a size 3 or 4:

Sizes - Circumferences in cm			
2 (II)	3 (III)	4 (IV)	5 (V)
30-37	33-40	35-43	37-46
20-22	22-24	24-26	26-28

If the measurements fall between two sizes use your clinical judgement for each patient.



Length - mediven for men	
ID regular	39-44
ID petite	33-38

Length - mediven for men	
ID regular	43-49
ID petite	37-42

Sizes							
Circumferences in cm	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)
cC	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cB	18-20	20-22	22-24	24-26	26-28	28-30	30-32

PIP CODES:

Please refer to our brochures or website (www.mediuk.co.uk) for PIP Codes.