


# mediven® Ulcer Kit

## Patient Referral Form FP10/GP10

mediven Ulcer Kit is listed under Part IXA of the Drug Tariff, under 'Venous Ulcer Compression System'

From: _____	
Patient's name: _____	
Clinician: _____	Date: _____
Contact Tel: _____	

Your patient should present this sheet to their GP to obtain product on Prescription.



**IA-D** Length of garment following contour of leg

**TO MEASURE:**

- Select length required. Measure A (floor) to D (below knee)  
>38.5cm = standard, <38.5cm = petite.
- Measure circumference at B (ankle at narrowest point), C (calf at widest point) and D (below knee).
- Select size required and identify appropriate code from table.

Sizes							
Circumferences in cm	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)
<b>cD</b>	27-33	29-36	32-39	34-42	36-45	38-48	40-50
<b>cC</b>	28-34	30-37	33-40	35-43	37-46	39-49	41-51
<b>cB</b>	18-20	20-22	22-24	24-26	26-28	28-30	30-32

mediven® Ulcer Kit Codes								
Product Code & Length		1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)
<b>Standard length &gt;38.5cm</b>	PIP Code	329-2109	329-2117	329-2125	329-2133	329-2141	329-2158	329-2174
	Product Code	DT275/I	DT275/II	DT275/III	DT275/IV	DT275/V	DT275/VI	DT275/VII
	<b>Quantity</b>							
<b>Petite length &lt;38.5cm</b>	PIP Code	329-2190	329-2208	329-2216	329-2224	329-2232	329-2240	329-2257
	Product Code	DT276/I	DT276/II	DT276/III	DT276/IV	DT276/V	DT276/VI	DT276/VII
	<b>Quantity</b>							
<b>Liner Socking</b>	PIP Code	329-2265	329-2273	329-2281	329-2299	329-2307	329-2315	329-2323
	Product Code	DT277/I	DT277/II	DT277/III	DT277/IV	DT277/V	DT277/VI	DT277/VII
	<b>Quantity</b>							

The mediven Ulcer Kit pack contains: 2 liner stockings, 1 over stocking.  
Available in 2 lengths and 7 sizes.  
Spare liner stockings are available in packs of 2 (1 length & 7 sizes).

For further information, please contact: medi UK Ltd, Plough Lane, Hereford, HR4 0EL.  
Tel: 01432 373500 | enquiries@mediuk.co.uk | www.mediuk.co.uk

