



mediven® harmony armsleeve & glove (RAL ccl 1 & 2)

RAL standard round knit armsleeves and gloves available on prescription

This form is designed to help you to identify which harmony product suits your patient needs and to ensure that all details are available for correct prescribing.

From:

Patient's name:

Measured by: Date:

Contact tel. no.

Please tick relevant boxes to identify product required

Armsleeve colour: Standard = Sand | Short = Caramel
 Glove colour: Sand

1 mediven® harmony

2 RAL Compression Class: ccl 1 ccl 2

3 Length: standard short

4 Style:

<input type="checkbox"/>	Armsleeve (C-G)	<input type="checkbox"/>	Gloves (sizes 2 - 7)
<input type="checkbox"/>	Armsleeve Extra Wide (C-G)		
<input type="checkbox"/>	** Armsleeve with Grip Top (C-G)		
<input type="checkbox"/>	** Armsleeve with Extra Wide Grip Top (C-G)		
<input type="checkbox"/>	Armsleeve with Shoulder Cap (C-G)		
<input type="checkbox"/>	Armsleeve Extra Wide with Shoulder Cap (C-G)		
<input type="checkbox"/>	Combined Armsleeve (A-G)		
<input type="checkbox"/>	Combined Armsleeve Extra Wide (A-G)		
<input type="checkbox"/>	** Combined Armsleeve with Grip Top (A-G)		
<input type="checkbox"/>	** Combined Armsleeve Extra Wide with Grip Top (A-G)		
<input type="checkbox"/>	Combined Armsleeve with Shoulder Cap (A-G)		
<input type="checkbox"/>	Combined Armsleeve Extra Wide with Shoulder Cap (A-G)		

**** Short armsleeve length available in these styles**

5 Size: 1 (I) 2 (II) 3 (III) 4 (IV) 5 (V) 6 (VI) 7 (VII) 8 (VIII)


Product Codes (if known):	PIP Codes (if known):	Quantity Required:
---------------------------	-----------------------	--------------------

mediven® harmony measurement chart (RAL ccl 1 & 2)



Armsleeve measurements

Circ (cms)	Size 1 (I)		Size 2 (II)		Size 3 (III)		Size 4 (IV)		Size 5 (V)		Size 6 (VI)		Size 7 (VII)		Size 8 (VIII)	
	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide
cG	25.5 - 28	28.5 - 31	28 - 30.5	30.5 - 33.5	30 - 33	33 - 36	32 - 35	35 - 38.5	33.5 - 36.5	37 - 40.5	36 - 39.5	39.5 - 43	38 - 41.5	41.5 - 45.5	40.5 - 44.5	45 - 49.5
cE	23 - 25	24 - 26	24.5 - 26.5	26.5 - 29	26 - 28.5	28.5 - 31	27.5 - 30	30 - 33	29 - 31.5	31.5 - 34.5	30.5 - 33.5	33.5 - 36.5	32.5 - 35.5	35.5 - 39	35 - 38.5	38.5 - 42
cC	13.5 - 14.5		14.5 - 15.5		16 - 17.5		17.5 - 19		19 - 20.5		20.5 - 22.5		22.5 - 24.5		25 - 27.5	
cA	15 - 18		16 - 19		17 - 20		19 - 22		20 - 23		22 - 25		24 - 27		26 - 29	



Measuring points

Armsleeve	C, E & G
Combined Armsleeve	A, C, E & G
Glove	A & C

Armsleeve lengths (C-G) (cms)

Standard	41-47
Short	<41

Glove measurements

Circ (cms)	Size 2 (II)	Size 3 (III)	Size 4 (IV)	Size 5 (V)	Size 6 (VI)	Size 7 (VII)
cC	13.5 - 15.5	15 - 18	17 - 20	19 - 21	20 - 23.5	22.5 - 27.5
cA	15 - 19	17 - 20	19 - 22	20 - 23	22 - 25	24 - 29

medi help and advice line - T. 01432 373500

medi UK Ltd, Plough Lane, Hereford, HR4 0EL - T. 01432 373500
 enquiries@mediuk.co.uk - www.mediuk.co.uk - shop.mediuk.co.uk