

From: .....

Patient's name: .....

Measured by: ..... Date: .....

Contact Tel No: .....



**medi UK Ltd**  
Plough Lane  
Hereford  
HR4 0EL  
Tel. 01432 373500  
Fax. 01432 373510  
www.mediuk.co.uk


### Juxta-Fit Essentials Upper Leg with knee piece

Product & PIP Codes	Size circumference (cms)			35cm Length - Right		35cm Length - Left		45cm Length - Right		45cm Length - Left		55cm Length - Right		55cm Length - Left	
	cD	cE1	cG		√		√		√		√		√		√
<b>XS</b>	<b>27-35</b>	<b>37-48</b>	<b>46-62</b>	DT 954R/1 PIP: 404-1216		DT954L/1 PIP: 404-1133		DT 955R/1 PIP: 378-6548		DT 955L/1 PIP: 378-6498		DT 956R/1 PIP: 378-6654		DT 956L/1 PIP: 378-6597	
<b>S</b>	<b>31-41</b>	<b>42-55</b>	<b>54-70</b>	DT 954R/2 PIP: 404-1224		DT954L/2 PIP: 404-1117		DT 955R/2 PIP: 378-6555		DT 955L/2 PIP: 378-6506		DT 956R/2 PIP: 378-6662		DT 956L/2 PIP: 378-6605	
<b>M</b>	<b>35-46</b>	<b>50-65</b>	<b>60-80</b>	DT 954R/3 PIP: 404-1232		DT954L/3 PIP: 404-1158		DT 955R/3 PIP: 378-6563		DT 955L/3 PIP: 378-6514		DT 956R/3 PIP: 378-6670		DT 956L/3 PIP: 378-6613	
<b>L</b>	<b>40-52</b>	<b>58-76</b>	<b>69-90</b>	DT 954R/4 PIP: 404-1240		DT954L/4 PIP: 404-1174		DT 955R/4 PIP: 378-6571		DT 955L/4 PIP: 378-6522		DT 956R/4 PIP: 378-6688		DT 956L/4 PIP: 378-6621	
<b>XL</b>	<b>44-58</b>	<b>66-87</b>	<b>77-101</b>	DT 954R/5 PIP: 404-1257		DT954L/5 PIP: 404-1190		DT 955R/5 PIP: 378-6589		DT 955L/5 PIP: 378-6530		DT 956R/5 PIP: 378-6696		DT 956L/5 PIP: 378-6639	


### Juxta-Fit Lower Legging / Lower Legging Extra Wide Calf

Product & PIP Codes	Size circumference (cms)			Petite Length		Standard Length	
	cB	cB1	cC	28cm (11")	√	36cm (14")	√
<b>S</b>	15-23	20-28	30-38	DT 950/2 PIP: 348-0548		DT 951/2 PIP: 348-0597	
<b>M</b>	20-28	25-33	36-43	DT 950/3 PIP: 348-0555		DT 951/3 PIP: 348-0605	
<b>M-XW</b>	20-28	33-40	43-50	DT 950/3X PIP: 368/0782		DT 951/3X PIP: 368-0816	
<b>L</b>	25-33	30-38	41-48	DT 950/4 PIP: 348-0563		DT 951/4 PIP: 348-0613	
<b>L-XW</b>	25-33	38-46	48-56	DT 950/4X PIP: 368-0790		DT 951/4X PIP: 368-0824	
<b>XL</b>	30-38	36-46	46-56	DT 950/5 PIP: 348-0571		DT 951/5 PIP: 348-0621	
<b>XL-XW</b>	30-38	41-51	51-61	DT 950/5X PIP: 368-0808		DT 951/5X PIP: 368-0832	
<b>XXL</b>	33-41	43-53	53-64	DT 950/6 PIP: 348-0589		DT 951/6 PIP: 348-0639	


Upper leg measurements:



Lower leg measurements :



Foot measurements:



**Juxta-Fit Lower Leg** - Product length: petite 28cm or standard 36cm. Measured from medial ankle bone to knee crease in standing position.

**Juxta-Fit Ankle Foot Wrap** - Closed heel, foot length from base of toes 16cm / Inter-locking (open heel), foot length from base of toes 13-25cm.

### Juxta-Fit Ankle Foot Wrap, Interlocking (open) & Closed Heel

Product & PIP Codes	Size circumference (cms)			Open Heel		Closed Heel	
	length A-Y	cH	cB		√		√
<b>S</b>	>9	19-24	>36	DT 960/2 PIP: 348-0647		DT 961/2 PIP: 348-0670	
<b>M</b>	>9	24-29	>36	DT 960/3 PIP: 348-0654		DT 961/3 PIP: 348-0688	
<b>L</b>	>9	29-34	>41	DT 960/4 PIP: 348-0662		DT 961/4 PIP: 348-0696	