

Juxtalite™



medi UK Ltd
 Plough Lane
 Hereford
 HR4 0EL
 Tel. 01432 373500
 Fax. 01432 373510

Patient Referral Form FP10/GP10

Your patient should present this sheet to their GP to obtain product on Prescription.


From:

Patient's name:

Clinician : Date:

Contact Tel.:

Juxtalite is listed under Part IXA of the Drug Tariff, under 'Venous Ulcer Compression System'



This garment is available in 2 lengths and 8 sizes.

To select the correct garment, measure length and circumferences as follows:

LENGTH:

- Length of garment from B (ankle) to D (2 fingers width below crease of the knee) following contour of leg. Select best fit.
28cm = short, 33cm = long.

CIRCUMFERENCES:

- Circumference at ankle (B) and widest calf (C).

L: Length of garment following contour of leg

Sizes								
Circ.	Small	Medium	Medium X Wide	Large	Large X Wide	X Large	X Large X Wide	XX Large
cC (cms)	26-36	33-44	44-54	41-51	51-61	37-48	48-61	53-64
cB (cms)	19-22	22-27		25-32		30-37		35-42

Juxtalite Codes									
Product Code & length		Small	Medium	Medium XW	Large	Large XW	X Large	X Large XW	XX Large
Short (28cm)	PIP Code	394-0996	394-1002	394-1044	394-1010	394-1051	394-1028	394-1069	394-1036
	Product Code	DT301/2	DT301/3	DT301X/3	DT301/4	DT301X/4	DT301/5	DT301X/5	DT301/6
	Quantity								
Long (33cm)	PIP Code	394-0913	394-0921	394-0962	394-0939	394-0970	394-0947	394-0988	394-0954
	Product Code	DT300/2	DT300/3	DT300X/3	DT300/4	DT300X/4	DT300/5	DT300X/5	DT300/6
	Quantity								

The Juxtalite pack contains 1 x Juxtalite, 1 x BPS card, 2 x Comfort leg liners, 2 x Comfort compression anklets & 1 x Instructions for use.

For further information, please contact: medi UK Ltd, Plough Lane, Hereford, HR4 0EL.
 Tel: 01432 373500 | enquiries@mediuk.co.uk | www.mediuk.co.uk | shop.mediuk.co.uk