

physicool bandages comes in two sizes:

- Bandage Size A - Small (10cm x 2m) - suitable for wrists, elbows and ankles
- Bandage Size B - Large (12cm x 3m) - suitable for knees, thighs and shoulders
- The 500ml Coolant will recharge Bandage Size A and Bandage Size B.

PRODUCT	CODES	
	Single	Pack 25
Bandage Size A	PC01A	PC25A
Bandage Size B	PC01B	PC25B
	Single	Pack 15
Coolant 500ml	PC01C	PC15C



Bandage Size A



Bandage Size B



500ml Coolant



physicool 
Rapid Relief, Rapid Recovery

physicool is a unique cooling bandage,
combining cold, compression and support.

For more information or to purchase physicool products:

medi UK Ltd, Plough Lane, Hereford HR4 0EL.

Tel: 01432 373500 · E: enquiries@mediuk.co.uk

You can also purchase via our website: www.mediuk.co.uk





physicool is a CE marked cooling bandage that has been extensively used for First Aid in Sports Medicine since 2008. It produces rapid evaporation which works with the body to draw the heat away, rather than trying to drive cold in like most ice/cooling related products.

THE USE OF PHYSICOOOL FOLLOWING KNEE ARTHROPLASTY

physicool is currently being used within the NHS and private hospitals to reduce pain, swelling and inflammation across a wide range of applications.

Due to the initial success achieved after total knee arthroplasty; physicool is now being used to treat sprains, strains, bruising in tendons and ligaments making it an ideal trauma product.

physicool is also being used to treat joint inflammation, haematomas and rheumatoid arthritis, enhancing patient rehabilitation.

HOW PHYSICOOOL WORKS

physicool bandages come pre-treated with special Coolant. As the Coolant evaporates, heat is withdrawn from the limb.

The bandage is manufactured using an open-cell weave which acts as a conduit for the Coolant, it is 100% stretch cotton and contains no elastic. It is also finished with a Velcro strap that ensures that it cannot cause any pressure points.

Features and Benefits



- **INSTANT COOLING** - physicool's cooling properties work quickly, straight from the packet.
- **PRACTICAL AND MOBILE** - No need for ice or cold water. Patient free to mobilise while wearing. Patient can self apply.
- **PAIN REDUCTION** - Lowers the temperature of the skin and underlying issue, reducing patient pain.
- **SWELLING/INFLAMMATION REDUCTION** - Effective in reducing swelling and inflammation.
- **QUICKENS PATIENT MOBILITY** - Proven to help patients regain range of movement faster, speeding up mobility.
- **GOOD ANATOMICAL CONFORMITY** - When wrapped around an injured area, the bandage provides total coverage and all-round cooling and compression.
- **LONG LASTING** - Each bandage comes pre-treated and will provide between 2 to 3 hours of cooling before needing to be recharged using the 500ml Coolant.
- **REUSABLE** - The bandages are easily retreated, they are 100% cotton and washable.