

# Wrap Hunter<sup>®</sup> — to be used in conjunction with medi Hosiery Hunter<sup>®</sup>

medi

## Signs & Symptoms:

## Suitable juxta<sup>®</sup> wrap:

## Foot options:

All juxta products include compressive foot options

**Skin:** with or without a wound.  
**Limb shape:** good.  
**Oedema:** none to moderate (**juxtalite**), moderate to high (**juxtacures**).  
**Exudate levels:** none to moderate (**juxtalite**), moderate to high (**juxtacures**).

**juxtalite<sup>®</sup> pack**  
 1 x wrap | 2 x juxta undersocks  
 1 x BPS™ card | patient info leaflet  
**juxtacures<sup>®</sup> pack**  
 1 x wrap (body & spine) | 2 x comfort leg liners | 2 x compression anklets  
 1 x BPS™ card | Velcro security tabs  
 tape measure | patient info leaflet  
 \*(add in UCS™ to address skin & wound)

## Additional options:

### Are the toes swollen?

Yes = juxta toecaps



### Is more compression required to the dorsum of the foot?

Yes = interlocking foot wrap (open heel)



Yes = closed heel foot wrap for smaller feet / wheelchair users



## Venous disease

(patient assessed for suitable compression therapy)

## Chronic oedema / Lymphoedema

(patient assessed for suitable compression therapy)

### \*UCS™ Debridement

10 Sachets per box

Product Code DT500 | PIP Code 384-4271  
 NHSSC Code ELZ746 | PECOS Code 217106

**Skin:** normal or dry, hyperkeratotic, with or without wound.  
**Limb shape:** good or distorted with/without skin creases. Needs reduction or maintenance therapy.  
**Oedema:** mild to moderate.  
**Exudate levels:** mild to high including Lymphorrhoea (beading or trickling of fluid).

## juxtafit<sup>®</sup> lower legging pack


1 x wrap | 2 x juxta undersocks  
 1 x BPS™ card | 1 x silky cover-up  
 PAC Band | patient info leaflet  
 \*(add in UCS™ to address exudate & skin)

see reverse for size chart/order codes

see reverse for size chart/order codes

# Wrap Hunter<sup>®</sup> - measurement charts & codes:

## juxtalite:



Circumferences in cm	S (2)	M (3)	M - X (3X)	L (4)
cC	26 - 36	33 - 44	44 - 54	41 - 51
cB	19 - 22	22 - 27	22 - 27	25 - 32
	L - X (4X)	XL (5)	XL - X (5X)	XXL (6)
cC	51 - 61	37 - 48	48 - 61	53 - 64
cB	25 - 32	30 - 37	30 - 37	35 - 42
<b>Actual garment length in cm</b>	<b>short</b>	<b>long</b>		
lB - D	28	33		
<b>Product Codes</b>	<b>short</b>	<b>long</b>		
regular profile	DT 301/2-6	DT 300/2-6		
extra wide profile	DT 301X/3-5	DT 300X/3-5		

All garments contain the unique Built-In-Pressure system card (BPS).


## juxtacures:

### Circumferences in cm

Adjustable garments to cut according to individual circumferences.

	minimum	maximum		
cB	19	42		
cC	25	64		
<b>Actual garment</b>	<b>short</b>	<b>standard</b>	<b>long</b>	
lB - D	28	33	38	
<b>Product Codes</b>	<b>short</b>	<b>standard</b>	<b>long</b>	
	DT 100/2	DT 100/3	DT 100/4	
	<b>Comfort Liners</b>			
	DT 100L			

## juxtafit:





Circumferences in cm	S (2)	M (3)	M - X (3X)	L (4)
cC	30 - 38	36 - 43	43 - 50	41 - 48
cB1	20 - 28	25 - 33	33 - 40	30 - 38
cB	15 - 23	20 - 28	20 - 28	25 - 33
	L - X (4X)	XL (5)	XL - X (5X)	XXL (6)
cC	48 - 56	41 - 51	51 - 61	53 - 64
cB1	38 - 46	33 - 41	41 - 51	43 - 53
cB	25 - 33	30 - 38	30 - 38	33 - 41
<b>Actual garment length in cm</b>	<b>petite</b>	<b>regular</b>		
lB - D	28	36		
<b>Product Codes</b>	<b>petite</b>	<b>regular</b>		
regular profile	DT 950/2-6	DT 951/2-6		
extra wide profile	DT 950/3-5X	DT 951/3-5X		

## Interlocking foot wraps & juxta toe cap:

juxtafit interlocking ankle foot wrap	S	M	L
lA - Y	>9	>9	>9
cB	<36	<36	<41
cA1	19 - 24	24 - 29	29 - 34
lA	13 - 25	13 - 25	13 - 25
<b>Product Codes</b>	<b>S</b>	<b>M</b>	<b>L</b>
	DT 960/2	DT 960/3	DT 960/4

juxtafit ankle foot wrap	S	M	L
lA - Y	>9	>9	>9
cB	<36	<36	<41
cA1	19 - 24	24 - 29	29 - 34
lA	>16	>16	>16
<b>Product Codes</b>	<b>S</b>	<b>M</b>	<b>L</b>
	DT 961/2	DT 961/3	DT 961/4



juxta toe cap	2	3	4	5	6
cA Foot	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30
<b>Product Codes</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>
	DT 772/2	DT 772/3	DT 772/4	DT 772/5	DT 772/6