

# duomed soft® British Standard Hosiery - Patient Record

For patients **without oedema**.

All sizes below held in stock for next working day delivery.

Patient name: \_\_\_\_\_

Date: \_\_\_\_\_

Doctor/Clinician: \_\_\_\_\_

Pharmacist: \_\_\_\_\_

## Clinic or Pharmacy Address:


(Please specify the address where the patient will collect their hosiery)

**medi**

medi UK Ltd  
Plough Lane  
Hereford  
HR4 0EL  
Tel. 01432 373500  
Fax. 01432 373510  
www.mediuk.co.uk  
enquiries@mediuk.co.uk

## Request for duomed soft® British Standard Hosiery on FP10/GP10

For use as **Compression Liners**, select below knee, ccl1, open or closed toe



SIZE (circumference in cms)	SMALL (2)	MEDIUM (3)	LARGE (4)	X LARGE (5)	XX LARGE (6)
<b>cG THIGH</b>	42—57	48—64	54—71	60—78	66—85
<b>cC CALF</b>	28—34	32—38	36—42	40—46	44—50
<b>cB ANKLE</b>	19—21	22—24	25—27	28—30	31—34

1 PATIENT MEASUREMENTS	LEFT LEG	RIGHT LEG
cG 5cm below groin		
cC Widest part of calf		
cB Just above ankle bone		

2  = Size requested

3  = Number of pairs you require

### 4 Choose your compression class:

- BS Class 1**  
Light compression
- BS Class 2**  
Medium compression
- BS Class 3**  
Strong compression (open toe only)

### 5 Choose your style:

- Below knee**
- Thigh length**  
(with silicone topband)

### 6 Choose your colour:

- Sand**
- Black**

### 7 Toe option:


- Open toe**  
(includes applicator)
- Closed toe**
- 2in1 Hosiery Applicator**  
PIP Code: 329-8932 | Product Code: DT904

Your patient should present this sheet to their GP to obtain product on prescription

## duomed soft® British Standard Hosiery Range



### COMPRESSION CLASS 1 (ccl) - 14-17mmHg




BLACK - PIP CODES						
STYLE	PRODUCT CODE	SMALL (2)	MEDIUM (3)	LARGE (4)	X LARGE (5)	XX LARGE (6)
<b>Compression Liner</b>		Please order below knee, ccl1, open or closed toe				
Below knee, open toe	DT511/B/	388-1802	388-1810	388-1828	388-1836	388-1844
Below knee, closed toe	DT512/B/	388-1851	388-1869	388-1877	388-1885	388-1893
Thigh length with silicone topband, open toe	DT611/B/	388-1901	388-1919	388-1927	388-1935	388-1943
Thigh length with silicone topband, closed toe	DT612/B/	388-1950	388-1968	388-1976	388-1984	388-1992

### SAND - PIP CODES

STYLE	PRODUCT CODE	SMALL (2)	MEDIUM (3)	LARGE (4)	X LARGE (5)	XX LARGE (6)
<b>Compression Liner</b>		Please order below knee, ccl1, open or closed toe				
Below knee, open toe	DT511/	367-4850	367-4868	367-4876	367-4884	367-4892
Below knee, closed toe	DT512/	367-4926	367-4934	367-4959	367-4967	367-4975
Thigh length with silicone topband, open toe	DT611/	367-4983	367-4991	367-5006	367-5014	367-5022
Thigh length with silicone topband, closed toe	DT612/	367-5030	367-5048	367-5055	367-5063	367-5071

### COMPRESSION CLASS 2 (ccl) - 18-24mmHg



BLACK - PIP CODES						
STYLE	PRODUCT CODE	SMALL (2)	MEDIUM (3)	LARGE (4)	X LARGE (5)	XX LARGE (6)
Below knee, open toe	DT521/B/	388-2008	388-2016	388-2024	388-2032	388-2040
Below knee, closed toe	DT522/B/	388-2057	388-2065	388-2073	388-2081	388-2099
Thigh length with silicone topband, open toe	DT621/B/	388-2107	388-2115	388-2123	388-2131	388-2149
Thigh length with silicone topband, closed toe	DT622/B/	388-2156	388-2164	388-2172	388-2180	388-2198

### SAND - PIP CODES

STYLE	PRODUCT CODE	SMALL (2)	MEDIUM (3)	LARGE (4)	X LARGE (5)	XX LARGE (6)
Below knee, open toe	DT521/	367-5089	367-5097	367-5105	367-5113	367-5121
Below knee, closed toe	DT522/	367-5139	367-5147	367-5154	367-5162	367-5170
Thigh length with silicone topband, open toe	DT621/	367-5188	367-5196	367-5204	367-5212	367-5220
Thigh length with silicone topband, closed toe	DT622/	367-5238	367-5246	367-5253	367-5261	367-5279

### COMPRESSION CLASS 3 (ccl) - 25-35mmHg



SAND - PIP CODES						
STYLE	PRODUCT CODE	SMALL (2)	MEDIUM (3)	LARGE (4)	X LARGE (5)	XX LARGE (6)
Below knee, open toe	DT532/	367-5287	367-5295	367-5303	367-5311	367-5329
Thigh length with silicone topband, open toe	DT632/	367-5337	367-5345	367-5352	367-5360	367-5378

For further information, please call 01432 373 500 to speak directly with our Customer Services Team or email [clinical.team@mediuk.co.uk](mailto:clinical.team@mediuk.co.uk)