


# mediven® harmony measurement chart (RAL ccl 1 & 2)



## Armsleeve measurements

Circ (cms)	Size 1 (I)		Size 2 (II)		Size 3 (III)		Size 4 (IV)		Size 5 (V)		Size 6 (VI)		Size 7 (VII)		Size 8 (VIII)	
	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide	Regular	Extra Wide
cG	25.5 - 28	28.5 - 31	28 - 30.5	30.5 - 33.5	30 - 33	33 - 36	32 - 35	35 - 38.5	33.5 - 36.5	37 - 40.5	36 - 39.5	39.5 - 43	38 - 41.5	41.5 - 45.5	40.5 - 44.5	45 - 49.5
cE	23 - 25	24 - 26	24.5 - 26.5	26.5 - 29	26 - 28.5	28.5 - 31	27.5 - 30	30 - 33	29 - 31.5	31.5 - 34.5	30.5 - 33.5	33.5 - 36.5	32.5 - 35.5	35.5 - 39	35 - 38.5	38.5 - 42
cC	13.5 - 14.5		14.5 - 15.5		16 - 17.5		17.5 - 19		19 - 20.5		20.5 - 22.5		22.5 - 24.5		25 - 27.5	
cA	15 - 18		16 - 19		17 - 20		19 - 22		20 - 23		22 - 25		24 - 27		26 - 29	



### Measuring points

Armsleeve	C, E & G
Combined Armsleeve	A, C, E & G
Glove	A & C

### Armsleeve lengths (cms)

Standard	41-47
Short	<41

### Glove measurements

Circ (cms)	Size 2 (II)	Size 3 (III)	Size 4 (IV)	Size 5 (V)	Size 6 (VI)	Size 7 (VII)
cC	13.5 - 15.5	15 - 18	17 - 20	19 - 21	20 - 23.5	22.5 - 27.5
cA	15 - 19	17 - 20	19 - 22	20 - 23	22 - 25	24 - 29

**medi help and advice line - T. 01432 373500**

medi UK Ltd, Plough Lane, Hereford, HR4 0EL - T. 01432 373500  
 enquiries@mediuk.co.uk - www.mediuk.co.uk - shop.mediuk.co.uk