



**medi UK Ltd**  
 Plough Lane  
 Hereford  
 HR4 0EL  
 Tel. 01432 373500  
 Fax. 01432 373510

# mediven<sup>®</sup> Ulcer Kit

## Patient Referral Form FP10/GP10

Your patient should present this sheet to their GP to obtain product on Prescription.


From: .....

Patient's name: .....

Clinician : ..... Date: .....

Contact Tel.: .....

mediven Ulcer Kit is listed under Part IXA of the Drug Tariff, under '**Venous Ulcer Compression System**'



### Pack contents:

mediven Ulcer Kit - 2 liner stockings, 1 over stocking.  
 Available in 2 lengths and 7 sizes.  
 Spare liner stockings are available in packs of 2 (1 length & 7 sizes).

### To measure:

- Select length required. Measure A (floor) to D (below knee) >38.5cm = **standard**, <38.5cm = **petite**.
- Measure circumference at B (ankle at narrowest point), C (calf at widest point) and D (below knee).
- Select size required and identify appropriate code from table.

	Sizes						
Circ.	1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)
cD (cms)	27-33	29-36	32-39	34-42	36-45	38-48	40-50
cC (cms)	28-34	30-37	33-40	35-43	37-46	39-49	41-51
cB (cms)	18-20	20-22	22-24	24-26	26-28	28-30	30-32

mediven <sup>®</sup> Ulcer Kit Codes								
CODES & LENGTH		1 (I)	2 (II)	3 (III)	4 (IV)	5 (V)	6 (VI)	7 (VII)
Standard length >38.5cm	Product Code	DT275/I	DT275/II	DT275/III	DT275/IV	DT275/V	DT275/VI	DT275/VII
	PIP Code	329-2109	329-2117	329-2125	329-2133	329-2141	329-2158	329-2174
	Quantity							
Petite length <38.5cm	Product Code	DT276/I	DT276/II	DT276/III	DT276/IV	DT276/V	DT276/VI	DT276/VII
	PIP Code	329-2190	329-2208	329-2216	329-2224	329-2232	329-2240	329-2257
	Quantity							
Liner Stockings	Product Code	DT277/I	DT277/II	DT277/III	DT277/IV	DT277/V	DT277/VI	DT277/VII
	PIP Code	329-2265	329-2273	329-2281	329-2299	329-2307	329-2315	329-2323
	Quantity							

For further information, please contact: medi UK Ltd, Plough Lane, Hereford, HR4 0EL.  
 Tel: 01432 373500 | enquiries@mediuk.co.uk | www.mediuk.co.uk | shop.mediuk.co.uk

