

A woman with her hair in a ponytail, wearing a white towel around her neck and a blue wrist brace on her right hand, is sitting on a gym machine and holding a clear plastic water bottle. The background shows a blurred gym environment with various pieces of equipment.

medi

Therapy in safe hands – the medi hand programm

medi. I feel better.





Manumed tri wrist orthosis to restrict pronation and supination

medi hand products provide the right medical care for each indication and offer a high level of technology for maximum customer satisfaction.

Your key medi benefits

- medi products are easy to handle, low profile and light weight
- with medi compression technology and other benefits, medi products provide premium technological standards
- medi is constantly researching and applying new technologies
- high compliance and customer satisfaction confirm medi as a leading global player in the healthcare market

Manumed tri


Wrist brace for limiting pronation and supination

Indications

- carpal tunnel syndrome
- tenovaginitis
- restricted freedom of movement after injuries of the triangular fibrocartilage complex
- condition following fracture of the radius for functional therapy
- restricted freedom of movement after disc injuries
- osteoarthritis of the distal radioulnar joint
- guyon's canal syndrome
- condition after wrist fractures
- condition after disc surgery / disc suture

Mode of action

- Limitation of movement of the wrist in three planes of movement:
- Flexion / extension
- Abduction / adduction
- Pronation / supination
- Limitation of the range of movement of the DRUG



Manumed tri

size	circumference wrist cm	article number	article number
		left	right
I	14 – 18	R.111.201	R.111.301
II	19 – 23	R.111.202	R.111.302



Limitation of pronation



Neutral position



Limitation of supination

Features



Open frame construction Makes it very

comfortable to wear, patients do not sweat and can be worn at night. The orthosis is simple to use both before and after the operation; this avoids double prescriptions.



Optimum fit

The aluminium construction enables it to be fitted simply and individually to the patient's hand.



Higher wearing comfort

Comfort zone in the thumb region aids fit and compliance.



More mobility for the patients in everyday living

Innovative frame construction enables the hand to be held in the gable position.



Simple handling

The slip-on system with integrated securing loops ensure easy application.

Manumed active

Soft wrist support with stabilising splint

Indications

- wrist injuries, e.g. sprains, strains
- osteoarthritis of the wrist
- instability
- tenosynovitis
- irritative conditions (post-traumatic, postoperative and rheumatic)

Mode of action

- the volar splint stabilises and reduces stress on the wrist
- the circular strap complements the stress relief provided by the compressive knit and the volar hand splint
- pinpoint pressure relief on the ulnar styloid process enhanced by the anatomically designed strap



Coloures



silver



sand



Manumed active

size	circumference wrist cm	style	article number	
			left	right
I	13 – 15	silver	K.412.201	K.413.201
II	15 – 17	silver	K.412.202	K.413.202
III	17 – 19	silver	K.412.203	K.413.203
IV	19 – 21	silver	K.412.204	K.413.204
I	13 – 15	sand	K.412.801	K.413.801
II	15 – 17	sand	K.412.802	K.413.802
III	17 – 19	sand	K.412.803	K.413.803
IV	19 – 21	sand	K.412.804	K.413.804

Features



Pressure relief on the outer wrist joint with the strap opening that aids correct positioning.



Comfortable with highly elastic knitted weave in the thumb area combined with the ultra-flat seam.



Soft edge on forearm for added comfort, ensures a pressure-free edge without curling or constricting.




Effective support and stability with the optimally sized hand bar. The shape and dimensions of the hand bar are optimised for both small and large hands.




Pull-on aid supplied to assist with application - particularly useful if patient has sweaty hands or pain / inflammation in the area.

Features



Clima Comfort



Clima Fresh



medi compression – intelligent and effective products medi compression is the brand mark that makes medi products stand out. More than 60 years of German experience in design and research have been applied to bring you quality garments that are guaranteed to provide accurate graduated compression.

Manumed


Wrist support with stabilising splint

Indications

- chronic, post-traumatic or post-operative soft tissue injury of the wrist
- radial nerve paralysis
- immobilisation in distortions
- conservative treatment of tendonitis and/or tenosynovitis
- rheumatoid arthritis

Mode of action

- immobilisation and stabilisation of the wrist with mouldable aluminium splint




Features



Perfect fit achieved with anatomically formed, individually mouldable aluminium splint



Easy to fit 'slip-on' system



Comfort guaranteed with medi airtex material, an innovative, breathable, light and compressive material for comfortable activation of proprioception.

Manumed

size	circumference wrist cm	style	article number	
			left	right
0	up to 14	grey	R.110.230	R.110.330
I	14-18	grey	R.110.231	R.110.331
II	19-23	grey	R.110.232	R.110.332
III	24-28	grey	R.110.233	R.110.333
0	up to 14	black	R.110.200	R.110.300
I	14-18	black	R.110.201	R.110.301
II	19-23	black	R.110.202	R.110.302
III	24-28	black	R.110.203	R.110.303

Manumed T


Wrist brace with thumb support

Indications

- sprain or distortion of the wrist and/or thumb
- postoperative fixation of the wrist and thumb
- lesion of the ulnar collateral ligament of the basal joint of the thumb
- tendinitis and tendinopathy of the wrist and thumb
- inflammatory phases of rheumatic diseases

Mode of action

- fixation and stabilisation of the wrist, the thumb saddle and the basal joint of the thumb.



Manumed T

size	circumference wrist cm	style	article number	
			left	right
0	up to 14	grey	R.110.260	R.110.360
I	14-18	grey	R.110.261	R.110.361
II	19-23	grey	R.110.262	R.110.362
III	24-28	grey	R.110.263	R.110.363
0	up to 14	black	R.110.220	R.110.320
I	14-18	black	R.110.221	R.110.321
II	19-23	black	R.110.222	R.110.322
III	24-28	black	R.110.223	R.110.323

Features




Increased mobility for the patient's everyday activities thanks to the innovative frame construction that permits flexion of the hand.



High wearing comfort with the Comfort Zone in the thumb area.



Easy to fit with the 'slip-on' system



Comfort guaranteed with medi airtex material, an innovative, breathable, light and compressive material for comfortable activation of proprioception.

Rhizomed®

Orthopaedic support for the thumb saddle joint

Indications

- injuries to the lateral ligaments of the thumb saddle joint (skier's thumb)
- irritations of the basic thumb joint and thumb saddle joint (e.g. rhizarthritis)
- distortions

Mode of action

- stabilisation of the thumb saddle joint and basic thumb joint



Features




High wearing comfort with an open construction that prevents both sweating and pressure on the surgical suture.



Single handed application with belt retaining loops.



Perfect fit with individually adjustable aluminium construction.



measurement: from the middle of the wrist to the end joint of the thumb (measured at the inner side of the hand)

Rhizomed

size	length cm	article number	article number
		left	right
S	8 – 10	R.140.202	R.140.302
M	10 – 12	R.140.203	R.140.303
L	over 12	R.140.204	R.140.304

medi CTS

Wrist support with finger splint

Indications

- conservative therapy and aftercare following surgery for carpal tunnel syndrome (CTS)
- immobilisation of wrist and fingers

Mode of action

- reduction of swelling, inflammation and pain by immobilisation of wrist and fingers




Features



safe and reliable immobilisation of the hand in a neutral position



Easy to fit 'slip-on' system.



measurement: from wrist to top of middle finger

medi CTS

size	cm	article number
XS	up to 15	R.110.000
S	15 – 17	R.110.001
M	17.5 – 19	R.110.002
L	19.5 – 22	R.110.003

medi wrist support

Wrist support

Indications

- chronic, post-traumatic or post-operative soft tissue injury of the wrist
- radial nerve paralysis
- immobilisation after mild to moderate sprains
- conservative treatment of tendinitis and/or tenosynovitis
- rheumatoid arthritis
- after removing plaster cast

Mode of action

- immobilisation and stabilisation of the wrist



medi wrist support			
size	circumference wrist cm	article number	article number
		left	right
S	10 – 14	R.200.201	R.200.301
M	14 – 18	R.200.202	R.200.302
L	18 – 22	R.200.203	R.200.303
XL	22 – 26	R.200.204	R.200.304

medi thumb support


Thumb support

Indications

- injuries of the carpometacarpal and metacarpophalangeal joints of the thumb
- ligamentous laxity of the thumb
- arthrosis of the carpometacarpal joint of the thumb
- rheumatoid arthritis

Mode of action

- immobilisation of the thumb



medi thumb support			
size	circumference wrist cm	article number	article number
		left	right
S	10–14	R.210.201	R.210.301
M	14–18	R.210.202	R.210.302
L	18–22	R.210.203	R.210.303
XL	22–26	R.210.204	R.210.304

Headquarters
medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth
Germany
T + 49 921 912-0
F + 49 921 912 781
export@medi.de
www.medi.de



9E002/04.2013