



medi UK Ltd
 Plough Lane
 Hereford
 HR4 0EL
 Tel. 01432 373500
 Fax. 01432 373510

mediven[®] RAL hosiery

- mediven[®] for men
- mediven[®] active

Selecting the correct garment for your patient.

Your patient should present this sheet to their GP to obtain product on Prescription.

From:

Patient's name:

Measured by: Date:

Contact Tel.:

mediven[®] for men

mediven[®] active

Closed Toe

Closed Toe

Choose your RAL Compression Class:

Class 1 Class 2

18-21mmHg 23-32mmHg

Choose your RAL Compression Class:

Class 1 Class 2

18-21mmHg 23-32mmHg

Choose your style:

Length: Standard Petite

Below Knee (pair)

Standard length ID (see reverse) - 39-44cm
 Petite length ID (see reverse) - 34-49cm

Order Codes - Class 1		Order Codes - Class 2	
Standard	Petite	Standard	Petite
DT342	DT343	DT340	DT341

Choose your style:

Length: Standard Petite

Below Knee (pair)

Standard length ID (see reverse) - 43-49cm
 Petite length ID (see reverse) - 37-42cm (black only)

Order Codes - Class 1		Order Codes - Class 2	
Standard	Petite	Standard	Petite
DT180	DT181	DT280	DT281

Choose your colour:

Black (.../5)
 Navy (.../3)
 Grey (.../2)

Choose your colour:

Black (.../5) White (.../7)
 Navy (.../3) Brown (.../6)
 Grey (.../2)

Choose your size:

I II III IV V VI VII


Pairs Required:

DT Code if known:

See over for measurement chart



MEASURING:



calf length

A-D = length cD = circumference

Length mediven for men

∠ D regular	39 – 44
∠ D petite	34 – 38

Length mediven active

∠ D regular	43 – 49
∠ D petite	37 – 42

MEASUREMENT CHART:

Circumference (cm)

Sizes	I	II	III	IV	V	VI	VII
calf	eC 28 – 34	30 – 37	33 – 40	35 – 43	37 – 46	39 – 49	41 – 51
ankle	eB 18 – 20	20 – 22	22 – 24	24 – 26	26 – 28	28 – 30	30 – 32

Example:

Circumference B = 24cm, circumference C = 40cm, therefore could be a size 3 or 4:

Circumference: size in cm

Size II III IV V

eC 30 – 37 33 – 40 35 – 43 37 – 46

eB 20 – 22 22 – 24 24 – 26 26 – 28

- 1) Measure the circumference at the narrowest part of the ankle (B).
- 2) Measure the circumference at the widest part of the calf (C).
- 3) Measure the length of the limb at the back of the heel (flat on the floor) to a point 2 finger widths down from the knee flexure.

ORDER CODES & PIP CODES:

Please refer to Formulary Guide or website (www.mediuk.co.uk) for Order & PIP Codes.

If the measurements fall between two sizes use your clinical judgement for each patient.